


MICRODOME



Bipolar Coagulation Electrode



FOCUS ON BETTER SOLUTIONS




- 01 Sehr feiner bipolarer Elektrodenkopf für präzise Koagulationen.
Very fine bipolar electrode tip for precise coagulations.
- 02 Biegbare Version für mehr Flexibilität während der Anwendung.
Bendable version for more flexibility during use.

DE Für Fälle in denen bipolare Pinzetten nicht mehr eingesetzt werden können, etwa zur bipolaren Koagulation diffus blutender Areale z.B. in der Nase, wurde deshalb eine neue Generation von Elektroden entwickelt, die mit einem Schaftdurchmesser von nur 2 mm dem Chirurgen ermöglichen auch in engsten Kavitäten eine effektive und sichere Blutstillung zu erreichen.

EN In cases where bipolar forceps cannot be used e.g. for bipolar coagulation of diffusely bleeding areas in the nose, a new generation of electrodes have been developed with a shaft diameter of only 2 mm. Microdome electrodes enable the surgeon to cauterise bleeding in an effective and safe way even in the smallest cavities.

Bipolare Koagulationselektrode MICRODOME
Bipolar Coagulation Electrode MICRODOME



Länge Length	Anschluss Plug	biegbar malleable	bajonett bayonet	gewinkelt 20° angled 20°
13 cm			85620410	85620420
14 cm		85620430		85620432
18 cm		85620460		85620462
24 cm		85620450		

max
250
Vp

137

95

1



Günter Bissinger
Medizintechnik GmbH
Hans-Theisen-Str.1
79331 Teningen/Germany

T +49 7641 9 14 33 0
F +49 7641 9 14 33 33
info@bissinger.com
www.bissinger.com

Bissinger – bessere Lösungen im Blick

Bipolare und monopolare Instrumente
Zubehör für die Elektrochirurgie
Chirurgische Instrumente
Entwicklung von Instrumenten auf
Kundenwunsch
Service und Reparaturen

Bissinger – Focus on better solutions

Bipolar and monopolar instruments
Accessories for electro-surgery
Surgical instruments
Design and development of instruments
to customer specifications
Service and repair