

89102002

Adapter bipolar



Technical specifications

Dimensions and Weights

Basic Specifications

Cable connector instrument side	Martin, Berchtold
Cable plug generator side	Valleylab, Lamidey, EMC

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 1 kVp
Packaging unit	1
GTIN 14	04250418736725