

89102005

Adapter monopolar



Technical specifications

Dimensions and Weights

Basic Specifications

Cable connector instrument side	Valleylab, Conmed, International
Cable plug generator side	Erbe, EMC/Dolley

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 5 kVp
Packaging unit	1
GTIN 14	04250418736756