

**89102005**

Adapter monopolar



## Technical specifications

### Dimensions and Weights

#### Basic Specifications

Cable connector instrument side	Valleylab, Conmed, International
Cable plug generator side	Erbe, EMC/Dolley

#### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

#### Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	N/A

#### Additional Product Information

Accessories rated voltage	max. 5 kVp
Packaging unit	1
GTIN 14	04250418736756