


85630201

Bipolar Electrode for MIS



Technical specifications

Dimensions and Weights

Total length	300 mm
Total length in cm	30 cm
Diameter	3 mm
Angle	30 °

Basic Specifications

Electrode shape	L-Hook
-----------------	--------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418739924