


85500003

## CLASSIC Bipolar Coagulation Forceps, Exterior Tube



### Technical specifications

#### Dimensions and Weights

Total length	340 mm
Diameter	5 mm

#### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

#### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

#### Additional Product Information

Packaging unit	1
GTIN 14	04250418722797