


85500007

CLASSIC Bipolar Coagulation Forceps, Exterior Tube



Technical specifications

Dimensions and Weights

Total length	450 mm
Diameter	5 mm

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Packaging unit	1
GTIN 14	04250418722827