


85500068

## CLASSIC Bipolar Coagulation Forceps, Electrode



### Technical specifications

#### Dimensions and Weights

Total length	200 mm
Diameter	5 mm

#### Basic Specifications

Electrode shape	Scissor with Tungsten Carbide Blade, curved
-----------------	---

#### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

#### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

#### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418723121