


**82710100**

## POWEREDGE Bipolar Coagulation Forceps, Electrode



## Technical specifications

### Dimensions and Weights

|              |        |
|--------------|--------|
| Total length | 340 mm |
| Diameter     | 5 mm   |

### Basic Specifications

#### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | Class Ir  |
| FDA Approval             | N/A       |

### Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 400 Vp    |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418737500 |