


82710104

POWEREDGE Bipolar Coagulation Forceps, Electrode



Technical specifications

Dimensions and Weights

Total length	250 mm
Diameter	5 mm

Basic Specifications

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	Class Ir
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 400 Vp
Packaging unit	1
GTIN 14	04250418740876