


**81517401**

**Bipolar Forceps, curved**



## Technical specifications

### Dimensions and Weights

Total length	170 mm
Total length in inches	6,69 "
Tip width	1 mm

### Basic Specifications

Forceps shape	curved
Tip shape	straight
Connector	Flat Plug

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418709033