


81122305

Bipolar Forceps Duckbill with Irrigation



Technical specifications

Dimensions and Weights

Total length	220 mm
Total length in inches	8,66 "
Tip width	8 mm

Basic Specifications

Forceps shape	straight
Tip shape	Duckbill
Connector	Flat Plug
with rinse	Yes

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418707695