81522351

Bipolar Forceps with Double Irrigation and Luer-Lock Connector





Technical specifications

Dimensions and Weights

Total length	220 mm
Total length in inches	8,66 "
 Tip width	1,5 mm

Basic Specifications

Forceps shape	bayonet
Tip shape	straight
Connector	Flat Plug
with rinse	pdf.checkbox.1

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
- Washing machine	95 °C
Autoclave sterilization	max. 137 °C
Accessories rated voltage	max. 500 Vp

Standards

Classification related to 93/42 EEC	Class IIb
FDA Classification	Klasse II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418710435