


81120303

Bipolar Forceps with Irrigation, bayonet



Technical specifications

Dimensions and Weights

Total length	200 mm
Total length in inches	7,87 "
Tip width	1,5 mm

Basic Specifications

Forceps shape	bayonet
Tip shape	straight
Connector	2-Pin Plug
with rinse	Yes

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418707602