


81122303

## Bipolar Forceps with Irrigation, bayonet



## Technical specifications

### Dimensions and Weights

Total length	220 mm
Total length in inches	8,66 "
Tip width	1,5 mm

### Basic Specifications

Forceps shape	bayonet
Tip shape	straight
Connector	2-Pin Plug
with rinse	Yes

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418707671