


81018203

Bipolar Forceps with Irrigation, straight



Technical specifications

Dimensions and Weights

Total length	180 mm
Total length in inches	7,09 "
Tip width	1,5 mm

Basic Specifications

Forceps shape	straight
Tip shape	straight
Connector	Flat Plug
with rinse	Yes

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418735476