


**81018204**

Bipolar Forceps with Irrigation, straight



## Technical specifications

### Dimensions and Weights

Total length	180 mm
Total length in inches	7,09 "
Tip width	0,5 mm

### Basic Specifications

Forceps shape	straight
Tip shape	straight
Connector	Flat Plug
with rinse	Yes

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418707268