


**81118202**

Bipolar Forceps with Irrigation, straight



## Technical specifications

### Dimensions and Weights

|                        |        |
|------------------------|--------|
| Total length           | 180 mm |
| Total length in inches | 7,09 " |
| Tip width              | 2 mm   |

### Basic Specifications

|               |            |
|---------------|------------|
| Forceps shape | straight   |
| Tip shape     | straight   |
| Connector     | 2-Pin Plug |

### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | Class II  |

### Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 500 Vp    |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418707510 |