


85791450

Bipolar Clamp-Scissor



Technical specifications

Dimensions and Weights

Total length	145 mm
Total length in cm	14,5 cm

Basic Specifications

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 250 Vp
Packaging unit	1
GTIN 14	04250418724579