


80100013

Bipolar Cable



Technical specifications

Dimensions and Weights

| | |
|--------------|-----|
| Cable length | 5 m |
|--------------|-----|

Basic Specifications

| | |
|---------------------------------|-----------|
| Cable connector instrument side | Flat Plug |
| Cable plug generator side | 2-Banana |

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|----------|
| Directive 93/42 EEC | N/A |
| Regulation (EU) 2017/745 | Class I |
| FDA Approval | Class II |

Additional Product Information

| | |
|---------------------------|----------------|
| Accessories rated voltage | max. 1 kVp |
| Packaging unit | 1 |
| GTIN 14 | 04250418700139 |