


80100045

Bipolar Cable



Technical specifications

Dimensions and Weights

Cable length	5 m
--------------	-----

Basic Specifications

Cable connector instrument side	2-Pin Plug
Cable plug generator side	Martin, Berchtold, Aesulap GK 55, GK 60

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 1 kVp
Packaging unit	1
GTIN 14	04250418700429