



**80100080**

**Bipolar Cable**



## Technical specifications

### Dimensions and Weights

Cable length	3 m
--------------	-----

### Basic Specifications

Cable connector instrument side	Round Plug for BiTech, MITHRAS, WAVE, PLASMALOOP
Cable plug generator side	Erbe, EMC/Dolley

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 1 kVp
Packaging unit	1
GTIN 14	04250418700610