



**80100084**

## Bipolar Cable



## Technical specifications

### Dimensions and Weights

Cable length	3 m
--------------	-----

### Basic Specifications

Cable connector instrument side	Round Plug for BiTech, MITHRAS, WAVE, PLASMALoop
Cable plug generator side	Martin, Berchtold, Aesulap GK 55, GK 60

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 1 kVp
Packaging unit	1

