


85717151

## BiTech Bipolar Scissor



## Technical specifications

### Dimensions and Weights

Total length	175 mm
Total length in cm	17,5 cm

### Basic Specifications

Scissor blade	Micro
Scissor handle	curved

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418733342