


85721120

BiTech Bipolar Scissor



Technical specifications

Dimensions and Weights

Total length	210 mm
Total length in cm	21 cm

Basic Specifications

Scissor blade	Delicate Blade
Scissor handle	straight

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418724463