

85713115

## BlackLine Bipolar Scissor



## Technical specifications

### Dimensions and Weights

|                    |        |
|--------------------|--------|
| Total length       | 130 mm |
| Total length in cm | 13 cm  |

### Basic Specifications

|                |          |
|----------------|----------|
| Scissor blade  | Micro    |
| Scissor handle | straight |

### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | Class II  |

### Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 300 Vp    |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418748285 |