


81450137

CLARIS+ Non-Stick Bipolar Forceps, bayonet



Technical specifications

Dimensions and Weights

Total length	180 mm
Total length in inches	7,09 "
Tip width	1,5 mm

Basic Specifications

Forceps shape	bayonet
Tip shape	straight
Connector	Flat Plug

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418741255