


81460158

CLARIS+ Non-Stick Bipolar Forceps, bayonet



Technical specifications

Dimensions and Weights

| | |
|------------------------|--------|
| Total length | 235 mm |
| Total length in inches | 9,25 " |
| Tip width | 0,6 mm |

Basic Specifications

| | |
|---------------|------------|
| Forceps shape | bayonet |
| Tip shape | straight |
| Connector | 2-Pin Plug |

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|-----------|
| Directive 93/42 EEC | Class IIb |
| Regulation (EU) 2017/745 | N/A |
| FDA Approval | Class II |

Additional Product Information

| | |
|---------------------------|----------------|
| Accessories rated voltage | max. 500 Vp |
| Packaging unit | 1 |
| GTIN 14 | 04250418738682 |