


**81460161**

**CLARIS+ Non-Stick Bipolar Forceps, bayonet**



## Technical specifications

### Dimensions and Weights

Total length	235 mm
Total length in inches	9,25 "
Tip width	1,5 mm

### Basic Specifications

Forceps shape	bayonet
Tip shape	straight
Connector	2-Pin Plug

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418741323