

**81450125**

**CLARIS+ Non-Stick Bipolar Forceps, angled**



## Technical specifications

### Dimensions and Weights

Total length	225 mm
Total length in inches	8,86 "
Tip width	0,3 mm

### Basic Specifications

Forceps shape	straight
Tip shape	angled
Connector	Flat Plug

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418739696