

80522332

CLARIS Non-Stick Bipolar Forceps, bayonet, angled



Technical specifications

Dimensions and Weights

Total length	220 mm
Total length in inches	8,66 "
Tip width	2 mm

Basic Specifications

Forceps shape	bayonet
Tip shape	angled up
Connector	Flat Plug

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418732154