

80611100

CLARIS Non-Stick Bipolar Forceps, Jeweler



Technical specifications

Dimensions and Weights

Total length	110 mm
Total length in inches	4,33 "
Tip width	0,25 mm

Basic Specifications

Forceps shape	Jeweler
Tip shape	straight
Connector	2-Pin Plug

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418705974