


80624337

CLARIS Xclusive Non-Stick Bipolar Forceps, bayonet



Technical specifications

Dimensions and Weights

Total length	240 mm
Total length in inches	9,45 "
Tip width	0,3 mm

Basic Specifications

Forceps shape	bayonet
Tip shape	straight
Connector	2-Pin Plug

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 250 Vp
Packaging unit	1
GTIN 14	0425041875170