


89101053

## Cable for Disposable Neutral Electrodes



## Technical specifications

### Dimensions and Weights

Cable length	4 m
--------------	-----

### Basic Specifications

Cable connector instrument side	for Disposable Neutral Electrode
Cable plug generator side	Valleylab (Non REM), International

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 5 kVp
Packaging unit	1
GTIN 14	04250418728003