

86620435

Mangeshikar Uterine-Manipulator Advanced+, Tube with Insulating Ring



Technical specifications

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|-----------|
| Directive 93/42 EEC | Class IIa |
| Regulation (EU) 2017/745 | N/A |
| FDA Approval | N/A |

Additional Product Information

| | |
|----------------|----------------|
| Packaging unit | 1 |
| GTIN 14 | 04250418749589 |