


93105100

Micro Dissection Needle Electrode, sterile



STERILE EO

Technical specifications

Dimensions and Weights

Total length	60,5 mm
Dimensions	mm
Diameter	0,3 mm

Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Needle

Cleaning and Reprocessing

Reusable	No
Manual reprocessing	No

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	20
GTIN 14	04250418737418