


93105110

Micro Dissection Needle Electrode, sterile



Technical specifications

Dimensions and Weights

Total length	57,3 mm
Dimensions	mm
Diameter	0,3 mm

Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Needle

Cleaning and Reprocessing

Reusable	No
Manual reprocessing	No

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	20
GTIN 14	04250418737456