


**89505100**

## Micro Dissection Needle Electrode



## Technical specifications

### Dimensions and Weights

Total length	60,5 mm
Dimensions	mm
Diameter	0,3 mm

### Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Needle

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

### Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	1
GTIN 14	04250418729093