

89505303

## Micro Dissection Needle Electrode



## Technical specifications

### Dimensions and Weights

|              |         |
|--------------|---------|
| Total length | 268 mm  |
| Dimensions   | mm      |
| Diameter     | 0,25 mm |

### Basic Specifications

|                 |        |
|-----------------|--------|
| Shaft diameter  | 2,4 mm |
| Electrode shape | Needle |

### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | N/A       |

### Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 2 kVp     |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418741774 |