



85530090

Mini Trocar 3.0 Micro, Sealing Cap



Technical specifications

Dimensions and Weights

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|---------|
| Directive 93/42 EEC | Class I |
| Regulation (EU) 2017/745 | N/A |
| FDA Approval | N/A |

Additional Product Information

| | |
|----------------|----------------|
| Packaging unit | 5 |
| GTIN 14 | 04250418745420 |