

82460312

MITHRAS Bipolar Coagulation Forceps, Electrode



Technical specifications

Dimensions and Weights

|              |        |
|--------------|--------|
| Total length | 150 mm |
| Diameter     | 3 mm   |

Basic Specifications

|                 |                                   |
|-----------------|-----------------------------------|
| Electrode shape | Flamingo Non-Stick Ø 0,7 mm / 45° |
|-----------------|-----------------------------------|

Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | Class II  |

Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 300 Vp    |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418746182 |