

82460340

MITHRAS Bipolar Coagulation Forceps, Electrode Bayonet Shape



## Technical specifications

### Dimensions and Weights

Total length	150 mm
Diameter	3 mm

### Basic Specifications

Electrode shape	Flamingo $\varnothing$ 0,7 mm / 15°
-----------------	-------------------------------------

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418746151