


85902402

Monopolar Wire Loop Electrode



Technical specifications

Dimensions and Weights

Diameter	125 mm
----------	--------

Basic Specifications

Electrode shape	Wire Loop
-----------------	-----------

Cleaning and Reprocessing

Reusable	No
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	1
GTIN 14	04250418735971