


89503006

Monopolar Electrode, 2,4 mm



Technical specifications

Dimensions and Weights

Total length	65 mm
Dimensions	mm
Diameter	0,7 mm

Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Needle

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 4 kVp
Packaging unit	1
GTIN 14	04250418728836