


89503009

Monopolar Electrode, 2,4 mm



Technical specifications

Dimensions and Weights

Total length	67 mm
Dimensions	2,5x20 mm
Diameter	0 mm

Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Knife

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 4 kVp
Packaging unit	1
GTIN 14	04250418728850