


**89503020**

Monopolar Electrode, 2,4 mm



## Technical specifications

### Dimensions and Weights

Total length	50 mm
Dimensions	mm
Diameter	2 mm

### Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Ball

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

### Additional Product Information

Accessories rated voltage	max. 4 kVp
Packaging unit	1
GTIN 14	04250418728911