


89503184

Monopolar Electrode, 2,4 mm



Technical specifications

Dimensions and Weights

Total length	150 mm
Dimensions	mm
Diameter	25 mm

Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Loop

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418729048