


89500150

Monopolar Electrode, 4 mm



Technical specifications

Dimensions and Weights

| | |
|--------------|-----------|
| Total length | 127 mm |
| Dimensions | 1,6x20 mm |

Basic Specifications

| | |
|-----------------|--------|
| Shaft diameter | 4 mm |
| Electrode shape | Lancet |

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|-----------|
| Directive 93/42 EEC | Class IIb |
| Regulation (EU) 2017/745 | N/A |
| FDA Approval | N/A |

Additional Product Information

| | |
|---------------------------|----------------|
| Accessories rated voltage | max. 4 kVp |
| Packaging unit | 1 |
| GTIN 14 | 04250418728690 |