

85900000

Monopolar Electrodes for MIS, Handle inclusive Shaft



Technical specifications

Dimensions and Weights

| | |
|--------------|--------|
| Total length | 320 mm |
| Diameter | 5 mm |

Basic Specifications

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|-----------|
| Directive 93/42 EEC | Class IIb |
| Regulation (EU) 2017/745 | N/A |
| FDA Approval | N/A |

Additional Product Information

| | |
|---------------------------|----------------|
| Accessories rated voltage | max. 2 kVp |
| Packaging unit | 1 |
| GTIN 14 | 04250418724647 |