

85900100

Monopolar Electrodes for MIS



Technical specifications

Dimensions and Weights

Diameter	5 mm
----------	------

Basic Specifications

Electrode shape	J-Hook
-----------------	--------

Cleaning and Reprocessing

Reusable	Yes
----------	-----

Manual reprocessing	Yes, related to IFU
---------------------	---------------------

Washing machine	95 °C
-----------------	-------

Autoclave sterilization	max. 137 °C
-------------------------	-------------

Approvals

Directive 93/42 EEC	Class IIb
---------------------	-----------

Regulation (EU) 2017/745	N/A
--------------------------	-----

FDA Approval	N/A
--------------	-----

Additional Product Information

Accessories rated voltage	max. 2 kVp
---------------------------	------------

Packaging unit	1
----------------	---

GTIN 14	04250418724708
---------	----------------